Bedienungsanleitung / Operating Instructions / Manuale d'uso

wenaFlex® Premium M

mit Kabelhandsender und Notabsenkung
with wired remote control and emergency lowering function
con telecomando a fili e abbassamento d'emergenza

Symbolfoto · Foto simbolica · Illustration is not binding

ACHTUNG:

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ATTENTION:

This device can be used by children aged 8 and over or by persons with reduced physical, sensory or mental capabilities or persons with lack of experience and knowledge, as long as they are supervised or have received instructions on safe use of the device and have understood potential resulting danger.

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English Contents

1.	Safety instructions	p.	9
2.	Startup and operation	p.	10
3.	Operation with wired remote control		
3.1.	Starting up	p.	11
3.2.	Wired remote control keys	p.	11
3.3.	Operating the bed adjustment functions	p.	11
3.4.	Emergency lowering	p.	12
4.	Technical data		
4.1.	Technical data on the drive	p.	12
4.2.	Technical data on the wired remote control	p.	12
4.3.	Mains cut-off	p.	12
5.	Troubleshooting	p.	13
6.	Disposal	p.	13
7.	Guarantee	p.	13
	APPENDIX: Engine removal	p.	19
	Declaration of conformity	p.	21

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1. Safety instructions

Safety instructions must be implicitly observed:

- In order to ensure safe and correct use, please read the operating instructions carefully before using the bed base.
- If the bed base is used incorrectly, it may present a potential danger to persons or cause damage to nearby objects.
- The bed base is only allowed to be used in dry rooms and has to be protected against exposure to moisture.
- The electrically adjustable bed base should not be kept outdoors.
- This electrically adjustable bed base is not licensed for use in the medical field.
- Electrical parts may only be opened or exchanged by qualified electricians.
- No persons or objects should be located beneath the raised parts of the bed base.
- **Electrically adjustable bed bases are not to be played with or used as toys by children.**
- **Small children should not be left unattended with electrically adjustable bed bases.**
- Cleaning and user maintenance tasks are not allowed to be performed by children without supervision.
- The drive system is only allowed to be used in continuous operation for a maximum duration of 2 minutes. Subsequent to this, the system needs a break of at least 18 minutes. If the system is used continuously for a duration of 1 minute, the system subsequently needs a break of at least 9 minutes.
- Before decommissioning or cleaning, disconnect the device from the mains power supply in order to eliminate the risk of electric shocks.
- Only use spare parts that are recommended by the manufacturer. Only such parts can ensure a sufficient level of safety.
- Used batteries are to be disposed of in an appropriate manner.

ATTENTION: This device can be used by children aged 8 and over or by persons with reduced physical, sensory or mental capabilities or persons with lack of experience and knowledge, as long as they are supervised or have received instructions on safe use of the device and have understood potential resulting danger.

2. Startup and operation

Please check and observe the following points:

- After storage or transport in cold conditions or during extreme temperature fluctuations, the device needs to be slowly adapted to the room temperature of the operating environment in order to avoid a development of water condensation before performing startup. If water condensation has developed, the device is not allowed to be put into operation for at least 12 hours.
- The transport safety devices are to be removed (i. e.: foil, carton, adhesive tape etc.)
- Are there any visible signs of transport damage on the bed base or on the drive/wired remote control? If so, please contact the Wenatex customer service department.
- Before performing the initial startup, check to ensure that the drive is securely fixed to the rod system. If not, please contact the Wenatex Customer Service Department.
- Is the bed construction suitable for this bed base? (See illustration 1)


Illustration 1

- Insert the bed base so that the frame has an even distance all the way around within the bed frame.
- Is the type of current and supply voltage in your home in compliance with the relevant technical data?
- Insert the plug into the electric socket.
- Check to ensure that the wired remote control is connected.
- Check to ensure that the mains connection and the control cable to the wired remote control are correctly wired.
- Protect all cables against crushing or other forms of potential damage. If any cables are crushed or damaged, immediately disconnect the mains plug from the power supply system!

- The electrically adjustable bed base is only to be used for the purpose of changing or adjusting the position of a person in a bed (see illustration 2).


Illustration 2

- While in a raised position, the head and foot sections are not suitable for sitting or standing on (see illustrations 3 and 4).

Illustration 3

Illustration 4

- Ensure that no persons or body extremities are located beneath or between the mechanical components of the bed base when raising or lowering any sections of the bed base. **Crushing hazard!**
- The electrically adjustable bed base is only intended for use as aforementioned. Any other form of use is inappropriate / inadmissible and can lead to accidents or the destruction of the device. Such inappropriate use will result in an immediate termination of any guarantee and warranty claims.
- The electrically adjustable bed base is not allowed to be used for lifting and lowering heavy objects in an industrial context or for work-related activities within the household.
- The electrically adjustable bed base is not allowed to be used in an environment in which it may be exposed to inflammable or explosive gasses or fumes (e.g. anaesthetics).
- The electrically adjustable bed base is not allowed to be used in damp environments (offshore installations, at or directly by the sea) or outdoors.
- The electrically adjustable bed base is not allowed to be used within aeroplanes or other forms of aircraft.
- The electrically adjustable bed base is not allowed to be used at nuclear power stations.

- Only use a dry or slightly damp cloth to clean the bed base (frame, motor, rod system, wired remote control).

Risk of danger to persons

- No liquids are allowed to enter the plug and socket connections.

Normal operation

- Observe the descriptions and the respective symbols of the wired remote control.
- Check the operation of the emergency lowering function every 6 months and then exchange the 9 volt battery. (Only use an alkaline battery!)

Battery-operated emergency lowering

If a power outage should occur while the adjustable bed base is in a raised position, you can lower it back down to a reclining position again with the aid of the wired remote control. The prerequisite for this is that the 9 volt battery is properly connected. **The capacity of the battery is only sufficient for one lowering sequence.** Exchange the battery after using for a lowering sequence.

3. Operation with wired remote control

3.1. Starting up

The system is delivered in a ready-to-use condition! Only the battery for the emergency lowering sequence needs to be connected to the system! Observe the following in order to avoid or remediate faults!

- Connect the battery for the emergency lowering sequence with the system. (This is absolutely necessary because the battery is also used for the mains cut-off.) (see illustration 7)
- Connect the wired remote control (A) to the socket (B) of the adjustment drive.
- Carefully pull the radio remote control cable through the strain relief clamp (C). (see illustrations 5 and 6).
- Finally, insert the plug with the mains cut-off (D) into the socket. (see illustration 5)
- The system is now ready for operation.

Illustration 5

Illustration 6

Illustration 7

3.2. Wired remote control keys

(see illustration 8)

- 1 Raise head section
- 2 Lower head section
- 3 Raise foot section
- 4 Lower foot section
- 5 Raise head and foot section simultaneously
- 6 Lower head and foot section simultaneously

Illustration 8

3.3. Operating the bed adjustment functions

Use the keys 1 to 6. When you press a key, the respective drive is set in motion. (The function LED indicator lights up when a key is pressed.)

3.4. Emergency lowering

The emergency lowering feature enables the system to retract to its basic position in case of a power failure. The drive is supplied with power via a 9V block battery. The battery is housed within the drive box. The battery is secured for transport and has to be connected to the battery clip before initial startup can be performed. (see information below point 3.1)

Note:

- Due to the severely limited capacity of the battery, only one single emergency lowering sequence can be performed. Then the battery has to be exchanged.
- Ensure that the battery is always fully charged.

Performance of emergency lowering

In order to perform emergency lowering, proceed as follows:

- Press the respective key(s) on the wired remote control that drive(s) the system back to the basic position.

4. Technical data

4.1. Twin drive incl. switch-mode power supply

Lifting force: **2 x 4,500 N**

Lift: **head 87 mm, foot 69 mm**

Motor voltage: **24 V DC motors**

Operating current: **3.15 A**

Power frequency: **50/60 Hz**

Power-on time: **10%, 2 min. on, 18 min. off**

Battery-operated emergency lowering:

9 V alkaline

Degree of protection: **IP20**

Protection class: **II**

Installation dimension: **581 mm, +/- 2 mm**

Handset socket: **5-pole**

4.2. Wired remote control

Input voltage: **29 V DC**

Control system: **Relay control**

Protection type: **IP20**

4.3. Mains cut-off

The mains cut-off feature prevents magnetic and electrical alternating fields on the bed, whereby no stand-by transformer is required.

Electrical current only flows when the hand control switch is activated, otherwise the bed remains de-energised.

A capacitor in the drive, which is charged with direct current, supplies power to the 2-pole relay in the mains cut-off and switches on the transformer in the drive. The capacitor is thereby recharged and is prepared for the next activation.

Note: When not in use over longer periods of time, the capacitor can lose its electrical potential. The installed 9 V battery of the emergency lowering feature then also buffers the capacitor of the mains cut-off.

Even in the unlikely case that the battery should also be empty, further operation can be performed without any problems. Press the green key (D) on the plug/mains cut-off (see illustration 5, page 11). The capacitor is then immediately charged and the mains cut-off is ready for operation again.

5. Troubleshooting – what to do when ...?

PROBLEM	CAUSE	SOLUTION
No voltage	Supply lines interrupted/ disconnected	Check supply line
No control voltage is provided	Fuses defect in the electrical distributor Customer's electrical socket defect	Exchange fuses Check/connect plug and socket connections
Wired remote control or drive system does not work	Wired remote control not plugged in	Check/connect plug and socket connections
	Drive system not plugged in	
	Wired remote control defect	Exchange wired remote control
Emergency lowering function cannot be performed	Monobloc battery empty	Exchange monobloc battery
	Objects are located beneath the bed base	Clear space
Suddenly the drive cannot be set in motion	Thermal circuit-breaker on the motor is possibly actuated (overload)	Leave the drive system in an inoperative state for approx. 20 to 30 minutes
	Capacitor for the mains cut-off is drained	Press the green key (D) on the plug/mains cut-off (illustration 5, page 11)


If the bed system still doesn't work properly, or if there appears to be damage to cables or connectors, please contact Wenatex Customer Service.

6. Disposal

The drive unit and the hand transmitter contain electronic components, cables, metals, plastics etc. The drive unit and the hand transmitter must be disposed in accordance with the applicable environmental regulations in your country.

The drive unit and the hand transmitter must not be disposed of with household waste. Batteries are not allowed to be disposed of along with household waste!

The instructions for removing the electrical components can be found on pages 19 and 20.

You can also find these operating instructions and safety instructions on:

www.wenatex.com

7. Guarantee

Wenatex provides a five year guarantee on all components of the Wenatex sleep system, starting as of the respective date of purchase. Fulfillment of a guarantee claim does not automatically suspend the running guarantee period, nor does it mark the beginning of a new guarantee period.

The following are excluded from the guarantee:

- Defects that were caused by inappropriate handling, i.e. use of force, overloading, changes performed by the buyer etc.
- Damages caused by opening up electrical or electronic components (unless the components were opened by a technician on behalf of the manufacturer).
- Defects caused by natural wear and tear, i.e. material fatigue of the spring elements, support slats or plastic parts.

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wenaFlex® Premium M bed base: Instructions for removing the engine

1. Disconnect the wenaFlex® Premium M bed base from the mains.
2. Remove the mattress to gain access to the wenaFlex® Premium M.
3. Disconnect the 9-volt battery for the emergency lowering from the drive unit.
4. Remove the marked spring elements (X) to gain access to the engine bracket (see Figure A). The spring elements are removed by turning them half way either clockwise or anticlockwise.
5. Slide both covers for the frame mount fully on and put them to one side (see Figure B).
6. Loosen the engine using pressure from above from the frame (see Figure C). The raised head and/or foot section should be held by a second person to reduce pressure from the frame. When removing the engine from the second frame, hold on to it to avoid damage to the floor and the engine.
7. **CAUTION: After removal of the engine, head and foot section fall back into the horizontal starting position (rigid position) (risk of jamming!)**

PLEASE NOTE: When head and foot section are raised to the maximum, the engine cannot be removed.




Please note the following dimensions: **Raising head section 425 mm (max. 575 mm); raising foot section 140 mm (max. 160 mm).**

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EC DECLARATION OF CONFORMITY

The manufacturer “Wenatex Forschung Entwicklung Produktion GmbH, Schloßstraße 49, 5282 Ranshofen” hereby declares that the following product: “electrically adjustable bed base”, product name: “wenaFlex® Premium M” meets all relevant provisions of the 2006/42/EC machinery directive.

The following harmonised standards were applied: **ÖNORM EN ISO 12100:2013**

Safety of machinery – General principles for design – Risk assessment and risk reduction

The following additional standards and guidelines were applied:

- 2014/30/EU electromagnetic compatibility
- 2014/35/EU low-voltage directive
- EN 60335-1:2012/A11:2014 safety of household and similar electrical appliances

The product is designed for the following operating conditions and operational environments:

Product may only be operated if it can be seen.

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Wenatex Das Schlafsystem GmbH
Wenatex Italia S.r.l.
Wenatex AG